


VE-1010

Veterinary PC ECG



EDAN

A world of potential

VE-1010 Veterinary PC ECG

EDAN's VE-1010 is a PC-based ECG solution for Veterinary, which includes a precise data sampling module and a suite of powerful PC software. It enables you to review, analyze, archive and network ECGs via your PC, aiming to provide you a paperless and streamlined ECG solution.

Features:

- Integrated data sampling box and software, paperless ECG solution for veterinaries
- Simultaneous 6/7 leads acquisition and recording
- Heart Rate Variability analysis
- Report printing on normal A4 paper
- Multi-formats reports: WORD, PDF, JPG and BMP
- Reports transmittable via email for teleconsultation possibility
- Networking your ECGs through LAN
- Complete report archives of patient information and ECG traces

