

GD 301 Evo

Constant heat medical sealer for the sealing of pouches and reels conforming to EN 868-5 and DIN 58953-7:2010 Directive.

Evolution of model GD 301, with improvements to the sealing unit and to the infeed and outfeed guides in order to increase performances: the sealer heats up much faster and more accurate, the film is perfectly kept flat inside the machine.

Reliable and user friendly, highly regarded for strong structure and attention to detail.

Compact machine, ideal for use in dentists, doctors and veterinary surgeries or laboratories.

Electronic control of sealing times to guarantee steady quality.

According to EN 868-5 / DIN 58953-7:2010

Sealing pouches up to 300 mm in width, manual actioning performed by a simple built in lever with an electronic timer giving consistent results every time.

The sealing system requires no consumables items, ensuring maintenance cost savings.

Vertically designed, including its integral roll holder, it can be positioned on any table even if there is limited width.

Accessories:

front support tray (included)

roll holder with cutting device (included)

MODEL	GD 301 Evo
Kind of sealing	Constant heat
Sealing width (mm)	12 flat
Max sealing length (mm)	300
Sealable materials	Pouches and reels EN 868-4
EN 8686-5 Conformity	Yes
DIN 58953-7:2010 Directive	Yes
Dimensions excl front support tray (mm)	452 x 286 x 242
Weight (kg)	3

